

**A FURTHER 5% OFF!**  
EXISTING OFFERS

# BLACK FRIDAY SALE



PROMOTION

NOV - DEC 2024



**MK-DENT**

ADAM  DENTAL  
supplies

# PRIME LINE



AeroAcoustic blade



Glass rod fiber optics



Reverse flow protection



Automatic pressure regulation<sup>1</sup>



Refined ceramic bearings



Scratch resistance  
Iconic Prime body

25  
W

Powerful & Precise  
Power Head - 25W

55  
dB

Incredibly Quiet  
Only 55 dB noise



Available for  
KaVo\*, NSK\*, W&H\*  
and Sirona\* Couplings



**BUY 1  
GET  
25% OFF**  
\$1,262.25<sup>EA</sup>

**BUY 2  
GET  
30% OFF**  
\$1,178.10<sup>EA</sup>

**BUY 3  
GET  
35% OFF**  
\$1,093.95<sup>EA</sup>

## Coupling Connection

Kavo  
NSK  
W&H  
Sirona

## Power Head 25 Watt

HP21KL  
HP21NL  
HP21WL  
HP21SL

## Small Head 18 Watt

HP22KL  
HP22NL  
HP22WL  
- N/A -

**\$1,683**



**2 YEARS  
UNCONDITIONAL  
WARRANTY\***

\*excluding physical damage on Prime Line handpieces

## COUPLINGS



- ✓ MK-dent & Kavo with LED & water
- QC6016KW

**\$495**



- ✓ NSK connection with LED
- QC6016NT

**\$385**



- ✓ W&H connection with fibre optic
- QC5016W

**\$385**



- ✓ Sirona connection with LED & water
- QC6016SW

**\$495**

- ✓ Non-optic
- QC4014K

**\$296**

- ✓ Non-optic
- QC4014NT

**\$296**

- ✓ Non-optic
- QC4014W

**\$296**

- ✓ Non-optic
- NA-

## ECO LINE TURBINES



- 62 dB Operating Noise
- AeroPerformance blade
- Refined ceramic bearings
- Glass rod fiber optics
- Smooth Eco body



Available for  
KaVo\*, NSK\*, W&H\*  
and Sirona\* Couplings



**BUY 2**  
**\$890**  
each

2 x EcoLine F/Optic



**BUY 3+**  
**\$792**  
each

3 x EcoLine F/Optic

Coupling Connection	Power Head 24W (F/Optic)	Small Head 16W (F/Optic)	Power Head 24W (N/Optic)	Small Head 16W (N/Optic)
KaVo	HE21KL	HE22KL	HE21K	HE22K
NSK	HE21NL	HE22NL	HE21N	HE22N
W&H	HE21WL	HE22WL	HE21W	HE22W
Sirona	HE21SL	- N/A -	- N/A -	- N/A -

**\$979 (FIBRE OPTIC)**

**\$693 (NON OPTIC)**



**1 Year FairCare Warranty**  
Unconditional warranty except physical damage

Buy 3  
**GET 10% OFF!**

## ECO LINE POWER EDITION



- 28 W Power
- 60 dB Operating Noise
- AeroPerformance blade
- Glass rod fiber optics
- Refined ceramic bearings
- Smooth Eco body

• HE20KL  
**\$1,265**

**BUY 3 or MORE**  
**\$1,089**  
each

Available for only  
KaVo Multiflex coupling fitting



**With 18 months FairCare Warranty**  
Unconditional warranty except physical damage

## ECO LINE FIXED CONNECTION TURBINES

- Direct 4 hole connection
- Refine ceramic bearings
- Push Button System
- HE17 - Standard head with 4 hole water spray
- HE18 - Small head with 3 hole water spray



**\$539**

Buy 3  
**GET 10% OFF!**

1) FairCare warranty (except physical damage) applicable on Eco Line sold after 1st Jan 2019 2) Eco Line Power Head 24 Watt

Buy 3  
**GET 20% OFF!**  
Mix & Match

**NEW  
RANGE!**



With 2 Years  
FairCare Warranty



SAFE-TECH  
PUSH  
BUTTON

⚙️ 1:5 Red ●

⌚ Max. 200,000 rpm

○ Internal water

⏏ For FG-burs  
(Ø 1.6 mm)

• LP15L (F/Optic)

**\$1,826**

**BUY 3**  
**\$1,460.80<sup>EA</sup>**



⚙️ 1:1 Blue ●

⌚ Max. 40,000 rpm

○ Internal water

⏏ For CA-burs  
(Ø 2.35 mm)

• LP11L (F/Optic)

**\$1,287**

**BUY 3**  
**\$1,029.60<sup>EA</sup>**



⚙️ 5.4:1 Green ●

⌚ Max. 7,407 rpm

○ Internal water

⏏ For CA-burs  
(Ø 2.35 mm)

• LP54L (F/Optic)

**\$1,320**

**BUY 3**  
**\$1,056<sup>EA</sup>**



⚙️ 1:1 Blue ●

⌚ Max. 40,000 rpm

○ Internal water

⏏ For HP-burs  
(Ø 2.35 mm)

• LP01L (F/Optic)

**\$1,045**

**BUY 3**  
**\$836<sup>EA</sup>**

## AIR MOTORS



**NEW**  
Shorter  
and  
Lighter

- ✓ Internal water
- ✓ LED

• AM1006

**\$990**



- ✓ Internal water

• AM1004

**\$748**



- ✓ External water

• AM0014

**\$473**



• AM0004

**\$473**

**NEW**

## KAVO COMPATIBLE HEADS



- Kavo 68LH compatible
- Blue (1:1)

• LCH11

**\$473**



- Kavo 67LH compatible
- Green (2:1)

• LCH21

**\$473**



- Kavo 31LR compatible
- Green (2:1)

• LCH21P

**\$484**

Buy 3  
**GET 1 FREE**  
\*of equal or lesser value  
**Mix & Match**



**FairCare Warranty<sup>1</sup>**  
1 year unconditional warranty  
except physical damage



SAFE-TECH  
PUSH  
BUTTON



⚙️ 1:5 Red ●  
○ Internal water  
For FG-burs  
(Ø 1.6 mm)

● LE15L (F/Optic)  
**\$1,595** **BUY 3 GET 1 FREE**

● LE15 (Non-optic)  
**\$1,375** **BUY 3 GET 1 FREE**

⚙️ 1:1 Blue ●  
○ Internal water  
For CA-burs  
(Ø 2.35 mm)

● LE11L (F/Optic)  
**\$1,155** **BUY 3 GET 1 FREE**

● LE11 (Non-optic)  
**\$891** **BUY 3 GET 1 FREE**

⚙️ 5.4:1 Green ●  
○ Internal water  
Screw-in & Snap-on  
Prophy Cups

● LE54LP (F/Optic)  
**\$1,287** **BUY 3 GET 1 FREE**

● LE54P (Non-optic)  
**\$1,045** **BUY 3 GET 1 FREE**

⚙️ 1:1 Blue ●  
○ Internal water  
For HP-burs  
(Ø 2.35 mm)

● LE01L (F/Optic)  
**\$990** **BUY 3 GET 1 FREE**

● LE01 (Non-optic)  
**\$825** **BUY 3 GET 1 FREE**

**Mix & Match**



⚙️ 1:1 Blue ●  
○ Push button  
● LB11/11

**\$418** **BUY 3 \$339<sup>EA</sup>**  
**BUY 4+ \$319<sup>EA</sup>**



⚙️ 1:1 Blue ●  
○ Latch Type  
● LB11/11C

**\$418** **BUY 3 \$339<sup>EA</sup>**  
**BUY 4+ \$319<sup>EA</sup>**

Buy More  
Save More



⚙️ 4:1 Green ●  
○ Snap-on Cups  
● LB41/11P

**\$594** **BUY 3 \$484<sup>EA</sup>**  
**BUY 4+ \$473<sup>EA</sup>**



⚙️ 1:1 Blue ●  
○ For HP-burs  
● LB02

**\$418** **BUY 3 \$339<sup>EA</sup>**  
**BUY 4+ \$319<sup>EA</sup>**

1) FairCare warranty (except physical damage) applicable on Eco Line sold after 1st Jan 2019



EMS Non-optic  
Type EN-041

- SC21E

\$462



EMS Piezon Non-optic  
Type EN-061

- SC21E2

\$561



EMS Piezon LED  
Type EN-060A

- SC21EL

\$748



KaVo Piezo LED  
Type 1.007.3995

- SC21KL\*

\$748



**NEW**



Acteon-Satelec Non-optic  
Type F 12 281

- SC21S

\$462



Acteon-Satelec LED  
Type F 12 609

- SC21SL

\$748



\*Only for use with EMS compatible scaler tips

## SCALER TIPS

**Buy 3 GET 1 FREE**  
While stocks last

- ✓ High-quality Diamond-Like-Carbon coating
- ✓ Extremely hard and wear-resistant surface

### • SCALING DLC COATED Only available for Satelec fitting



ST03SD  
General Scaling

**\$88**



ST04SD  
SATELEC 2

**\$132**

### • PERIO DLC COATED Only available for Satelec fitting



ST06SD  
SATELEC 10Z

**\$132**



ST10SD  
SATELEC 10P

**\$154**

### • ENDO DLC



ST08SD  
SATELEC S07

**\$154**

### • SCALING STAINLESS STEEL Only available for Satelec fitting



ST03S  
General Scaling

**\$63.80**



ST04S  
SATELEC 2

**\$88**

### • PERIO STAINLESS STEEL Available for EMS & Satelec fittings



ST07E / ST07S  
EMS DPL3  
SATELEC 10X

**\$110**



ST09E / ST09S  
EMS HPL

**\$110**



ST10E / ST10S  
SATELEC 10P

**\$110**



ST11E / ST11S  
EMS RT1

**\$110**

EMS and Acteon-Satelec are registered trademarks. MK-dent GmbH has no economic connection with the aforementioned companies.

# Universal Torque Wrench



TORQ CONTROL®

**Precision instrument with a Torque Adjustment System for tightening of Prosthetic components at specific Torque, for any implant system.**



## Automatic Declutching when Selected Torque is reached

- Compliance with Protocols for screw tightening
- Limits the risk of the screw fracturing or loosening



## Universal connection of Prosthetic mandrels

- 7 different Torque levels (from 10 to 35 N.cm)
- Can be used with ALL Implant Systems



## Micro-Head and Ergonomic Handle

- Easy access in restricted area with Excellent Visibility



## Easy and Quick to use

- Reduced motion compared with using a ratchet wrench while tightening



**\$1,210**


**Lightweight  
ONLY 135 gms**

Order code: 15501

## 3 Easy Steps Process



1. Please select the recommended Torque



2. Start rotating



3. Screw automatically stops when selected torque is reached



Made in France

[www.anthogyr.com](http://www.anthogyr.com)

**Anthogyr**  
A Straumann Group Brand

## Buy 1 Polisher

Get 1 Bottle of Propphy Powder FREE

Get 50% DISCOUNT ON 1 ADDITIONAL TIP

• **PR1011**

**\$1,375**

- Ergonomically shaped, compact and light-weight air-driven polishing device

- 360° rotatable outer handle due to unique RotaGrip® for simple and fatigue-proof operation



- Propphy Powder (300 gms) Strawberry - PP1013 \$60.50 each

**Additional tip**



• **SP3011**

**\$329**

**Coupling connection order codes:**

Kavo - PR1011K, NSK - PR1011N, W&H - PR1011W, Sirona - PR1011S

✓ **LEDs**

- ✓ Quality illumination, No Glare
- ✓ 5500K Natural Daylight Colour
- ✓ Clear Vision & Sharper Contrast
- ✓ Longer Durability & Economical



**\$88**

- BU8012 - KaVo & MK-dent coupling
- BU8012KM - KaVo & MK-dent motors
- BU8012WH - W&H handpieces
- BU8012WLT - W&H airmotors LT25
- BU8012WEA - W&H motors 40LT, 50 LT & A51LT
- BU8012SG - Sirona coupling and motors
- BU8012N - NSK coupling
- BU8012NT - NSK motors TiM40 & Ti205L
- BU8012NM - NSK coupling KCL,SCL,PTL LED

✓ **PREMIUM LUBRICATION OIL**



- ✓ Fully Synthetic
- ✓ For any brand of handpieces and motors

- LU1011 (500 ml)

**\$68.20**



- ✓ Fully Synthetic
- ✓ For KaVo Quattrocare 2 lubrication system

- LU1022 (500 ml)

**\$68.20**



- ✓ For W&H Assistina lubrication system

- LU1031 (500 ml)

**\$88**

MANUFACTURER



**MK-DENT**

MK-dent GmbH

Marie-Curie-Str. 2

22941 Bargteheide

GERMANY

[www.mk-dent.com](http://www.mk-dent.com)

YOUR SPECIALIST DEALER

**ADAM DENTAL**  
supplies

Ph: 1300 449 262

• All catalogue prices are inclusive of GST • Promotion valid until 31st Dec 2024 or while stocks last

KaVo, W&H, Sirona, NSK, Bien-Air, EMS and Acteon-Satelec are registered trademarks. MK-dent GmbH has no economic connection with the aforementioned companies.