



KULZER

MITSUI CHEMICALS GROUP

ADAM DENTAL
supplies

Giving a hand to oral health

Clinical Product Portfolio 2018/2019
Australia & New Zealand

Kulzer is your source for technologically advanced dental solutions that improve chair time productivity and effectiveness while revealing highly aesthetic results. Interrelated products like Venus Diamond®, Venus® Pearl, Venus White®, iBOND®, Gluma®, Flexitime® and Ivory® offer exceptional diagnosis, treatment, direct and indirect restorative care. This clinical product catalogue invites you to take a tour of our comprehensive portfolio of dental solutions.

www.kulzer.com

www.adamdental.com.au

Alginoplast	26
Charisma	9
Dispensing Impression Gun 1:1/2:1	22
Durafill VS	10
Dynamix Speed Unit	23
Fast & Scan / ProScan / Standard Flexitime	22-23
Flexitime BITE	25
Gluma Desensitizer	16
iBOND Systems	13
Impression Tips	22
Ivory System	20
Memosil 2 / Memoreg 2	25
Optosil / Xantopren	24
PLT Applicator	4
Sopira Citoject	17
Translux Wave light curing unit	18
Translux 2Wave light curing unit	19
Venus Bulk	8
Venus Colour	7
Venus Diamond	5
Venus Diamond Flow	6
Venus Pearl	4
Venus Supra Polishing Systems	17
Venus Whitening Systems	11
VibraKleen	27
Xantalgin	27
Xantasil	25

**Venus® Diamond and
Venus® Pearl**
Beauty beyond aesthetics.



RESTORATIVE SOLUTIONS

Venus® Pearl
Beauty beyond aesthetics.

Venus® Pearl is a universal nano-hybrid composite based on the innovative urethane monomer chemistry of Venus Diamond. Venus Pearl offers the same combination of low shrinkage stress and high durability—in a creamy consistency.

- Well suited for anterior as well as posterior restorations.
- Exceptional sculptability and high degree of polishability.
- Ideal for fine detailed work.

Innovative GUM & CORE shades with unique applications!



VENUS® PEARL

Venus Pearl PLT Master's Kit 10 x 0.2g PLT each of the shades A1, A2, HKA2.5, A3, A3.5, A4, 5 x 0.2g PLT each of the shades HKA5, B1, B2, B3, C1, C2, C3, D3, CL, AM, CO, YO, OLC, OMC, ODC, OXDC, BL, BXL, OB, GUM, 5 x 0.3g PLT of the shade CORE; 2 x 1.8g Venus Diamond Flow Syringe Baseline; Venus Family Shade Guide; Venus Pearl Layering Guide; Pictorial Card	KZ66048214
Venus Pearl Syringe Bundle Kit	KZ66065932
Venus Pearl PLT Bundle Kit	KZ66065931

ACCESSORIES	
1 PLT Applicator STEEL	KZ66019748

Shades	Syringe 3g	PLT 20 x 0.2g	PLT 10 x 0.2g
A1	KZ66048175	KZ66048090	
A2	KZ66048176	KZ66048141	
HKA2.5	KZ66048177	KZ66048142	
A3	KZ66048178	KZ66048143	
A3.5	KZ66048179	KZ66048144	
A4	KZ66048180	KZ66048145	
HKA5	KZ66048181	KZ66048146	
B1	KZ66048182	KZ66048147	
B2	KZ66048183	KZ66048148	
B3	KZ66048184	KZ66048149	
C1	KZ66048185	KZ66048150	
C2	KZ66048186	KZ66048161	
C3	KZ66048187	KZ66048162	
D3	KZ66048188	KZ66048163	

Shades	Syringe 3g	PLT 20 x 0.2g	PLT 10 x 0.2g
OLC (Opaque Light Chromatic)	KZ66048189		KZ66048164
OMC (Opaque Medium Chromatic)	KZ66048190		KZ66048165
ODC (Opaque Dark Chromatic)	KZ66048191		KZ66048093
BL (Bleach)	KZ66048193		KZ66048096
BXL (Bleach Extra Light)	KZ66048194		KZ66048097
OB (Opaque Bleach)	KZ66048195		KZ66048095
CL (Clear)	KZ66048196		KZ66048098
AM (Amber)	KZ66048197		KZ66048099
CO (Clear Opal)	KZ66048198		KZ66048100
YO (Yellow Opal)	KZ66048199		KZ66048171
CORE	KZ66048200		KZ66048172
GUM	KZ66048201		KZ66048173

RESTORATIVE SOLUTIONS

Venus® Diamond
Blending elegance and efficiency.

Venus® Diamond leads the industry in colour adaptability, with a patented formula that provides an outstanding natural look. Just one shade saves time and produces radiant results. Along with superb colour aesthetics, this revolutionary nano-hybrid composite exhibits exceptional durability, handling and minimal shrinkage.

- Revolutionary nano-hybrid that enables perfect colour adaptation, ease of polish and long lasting shine.
- Provides high durability with minimal shrinkage.
- Extended working time for efficient sculpting.

Use one shade or multiple shades to achieve a brilliant natural restoration!



VENUS® DIAMOND

Venus Diamond PLT Master's Kit (10 x 0.25g PLTs of the shades: A1, A2, HKA2.5, A3, A3.5, B1, B2, C2, OM, OD, 5 x 0.25g PLTs of the shades: A4, HKA5, B3, C3, D3, OL, CO, YO, CL, AM, BL, BXL, OB, 2 x 1.8g Venus Diamond Flow Baseline, and Shade Guide)	KZ66043143
--	------------

ACCESSORIES	
1 PLT Applicator STEEL	KZ66019748

Shades	Syringe 4g	PLT 20 x 0.25g	PLT 10 x 0.25g
A1	KZ66035471	KZ66039008	
A2	KZ66035472	KZ66039009	
HKA2.5	KZ66035473	KZ66039010	
A3	KZ66035474	KZ66039011	
A3.5	KZ66035475	KZ66039012	
A4	KZ66035476		KZ66039013
HKA5	KZ66035477		KZ66039014
B1	KZ66035478	KZ66039015	
B2	KZ66035479	KZ66039016	
B3	KZ66035480		KZ66039017
C2	KZ66035481	KZ66039018	
C3	KZ66035482		KZ66039019

Shades	Syringe 4g	PLT 20 x 0.25g	PLT 10 x 0.25g
D3	KZ66035483		KZ66039020
OLC (Opaque Light Chromatic)	KZ66046735		KZ66046731
OMC (Opaque Medium Chromatic)	KZ66046736	KZ66046732	
ODC (Opaque Dark Chromatic)	KZ66046737	KZ66046733	
BL (Bleach)	KZ66035487		KZ66039025
BXL (Bleach Extra Light)	KZ66035488		KZ66039026
OB (Opaque Bleach)	KZ66035489		KZ66039024
CL (Clear)	KZ66035490		KZ66039027
AM (Amber)	KZ66035491		KZ66039028
CO (Clear Opal)	KZ66035492		KZ66039029
YO (Yellow Opal)	KZ66035493		KZ66039030

RESTORATIVE SOLUTIONS

Venus® Diamond Flow

Seeing is believing.

Venus® Diamond Flow is a nano-technology based flowable composite that offers high radiopacity, low shrinkage and high flexural strength.

- Thixotropic handling behavior for easy control.
- Ideal as a baseliner.
- Unique Diamond chemistry produces superior esthetics through optimized color adaptation.



VENUS® DIAMOND FLOW

Venus Diamond Flow Syringe Assortment 1 x 1.8g syringe each of shades A1, A2, A3, and Baseliner	KZ66040353
--	------------

Venus Diamond Flow PLT Assortment 10 x 0.2g each of shades A1, A2, A3, and Baseliner	KZ66040367
---	------------

Shades	Syringe 1.8g	PLT 20 x 0.25g	Shades	Syringe 1.8g	PLT 20 x 0.25g
A1	KZ66040354	KZ66040368	HKA5	KZ66040361	KZ66040374
A2	KZ66040355	KZ66040369	B2	KZ66040362	KZ66040375
HKA2.5	KZ66040357	KZ66040370	OM (Opaque Medium)	KZ66040363	KZ66040376
A3	KZ66040358	KZ66040371	CL (Clear)	KZ66040364	KZ66040377
A3.5	KZ66040359	KZ66040372	BXL (Bleach Extra Light)	KZ66040365	KZ66040378
A4	KZ66040360	KZ66040373	Baseliner	KZ66041052	KZ66041051

RESTORATIVE SOLUTIONS

Venus® Color

As Individual as your Patients.

Many patients do not want their teeth to be recognised as new. They want their characteristic tooth-colour after treatment. Patients do not always want “brilliant white teeth”. A touch of colour is sometimes needed to restore your patient’s individual appearance.

We offer you Venus Colour in five different shades: White – Blue – Corn – Amber – Choco.
The drops of colour for aesthetics with an individual touch.

Whether recreating cracks in enamel, discoloured fissures or white spots, your restorations will look entirely natural. Venus Colour is the ideal accessory to our nano-hybrid composites Venus Diamond and Venus Pearl.

- Characterisation can be this simple.
- Venus Color is available in five shades, for natural-looking teeth.
- Easy and precise application from 1g syringes
- Suitable for direct and indirect fillings and prosthetics.
- Applicable in thin layers.

VENUS® COLOR

VENUS COLOUR SYRINGE ASSORTMENT Syringe Assortment	VENUS COLOUR SYRINGE REFILLS 1 x 1g	VENUS COLOUR ACCESSORIES Flow Cannulas
CONTENTS	Shades	Art. code
5 x 1g each of the shades WHITE. BLUE, CORN, AMBER, CHOCO 20 cannulas Venus Colour Shade Guide 25 brush tips Application instrument	WHITE	KZ66054206
	BLUE	KZ66054207
	CORN	KZ66054208
	AMBER	KZ66054209
	CHOCO	KZ66054210
Art. code	KZ66054221	




RESTORATIVE SOLUTIONS

Venus® Bulk Fill
Simply Beautiful. Simply Venus.

Venus® Bulk Fill is an advanced posterior restorative solution that provides dentists with a faster, easier technique than incremental layering. The result – a reduction in chair time, and a more efficient and cost effective restoration.

- Self-adaptive handling eliminates steps and provides ideal coverage of cavity floor and walls.
- Long working time for when you need it.
- Highly radiopaque for more accurate diagnosis.



VENUS® BULK FILL	
Venus Bulk Fill PLT Refill (10 x 0.20gm PLT's)	KZ66044743
Venus Bulk Fill Syringe Refill (1 x 1.8gm Syringe)	KZ66046795
Accessories	
1 PLT Applicator STEEL	KZ66019748

Venus Bulk Fill Specifications	
Curing Depth (blue light: 600 mWcm²; 20 sec)	6.2 mm
Shrinkage Stress (@24h)	3.4 MPa
Compressive Strength	331 MPa
Flexural Strength	120 MPa
Radiopacity	300%-Al
3 Body Wear Resistance: (ATCA method, poppy seed, 300,000 cycl.)	32 µm .31 mm
Working Time: (ambient light, 8klux)	80 seconds
Conversion Rate (blue light: 600 mWcm²; 20 sec)	55%
Reflection after Polishing with Venus Supra (20/40 sec)	7.5%

Source: Kulzer R&D, Wehrheim, Germany Unpublished data. Data on file

RESTORATIVE SOLUTIONS

Charisma®
Naturally Reliable.

The microparticle hybrid composite Charisma® has proven itself again and again over the last 24 years and is used every day by practitioners around the world. Millions of fillings have already been applied with Charisma worldwide. Thanks to Charisma's pronounced chameleon effect, patients' demands for aesthetically high-grade long-lasting restorations can be fulfilled. Now with ceramic reinforced formula.

- Charisma's high strength and wear resistance provide long-lasting restorations
- Charisma polishes to an extremely high surface luster
- Ultra-fine (0.7 microns) Microglass® particles result in a densely smooth surface comparable to natural enamel
- Available in 16 lifelike Vita®* shades for accurate shade matching and blending

CHARISMA®							
Shades	Syringe 4g	PLT 20 x 0.25g	PLT 10 x 0.25g	Shades	Syringe 4g	PLT 20 x 0.25g	PLT 10 x 0.25g
A1	KZ66000085	KZ66000103		C4	KZ66000095	KZ66000113	
A2	KZ66000086	KZ66000104		D3	KZ66000096	KZ66000114	
A3	KZ66000087	KZ66000105		OA2	KZ66000097		KZ66000115
A3.5	KZ66000088	KZ66000106		OA3	KZ66000098		KZ66000116
A4	KZ66000089	KZ66000107		OA3.5	KZ66000099		KZ66000117
B1	KZ66000090	KZ66000108		Accessories			
B2	KZ66000091	KZ66000109					
B3	KZ66000092	KZ66000110					
C2	KZ66000093	KZ66000111		Charisma Shade Guide		KZ66000121	
C3	KZ66000094	KZ66000112		1 PLT Applicator STEEL		KZ66019748	

*Vita® is a registered trademark of Vita Zahnfabrik GmbH



RESTORATIVE SOLUTIONS

Durafill® VS

The restorative choice for anterior restorations since 1979. Timeless Beauty.

Durafill® VS is a highly polishable and sculptable microfill resin for anterior restorations.

- Easy to place and contour—does not slump
- The standard to which other microfill composites are compared
- Excellent alternative for Class V restorations
- In addition to 16 Vita®* based shades, Durafill VS has Superlight shades, ideal for bleached, and primary teeth

DURAFILL® VS					
Shades	Syringe Refill	PLT Refill	Shades	Syringe Refill	PLT Refill
A1	KZ66000008	KZ66000027	OA2	KZ66000021	
A2	KZ66000009	KZ66000028	OA3	KZ66000022	
A3	KZ66000010	KZ66000029	OA3.5	KZ66000023	
A3.5	KZ66000011	KZ66000030	I	KZ66000026	
A4	KZ66000012		SL (Super light)	KZ66000025	
B1	KZ66000013		SSL (Super Snow Light)	KZ66002499	
B2	KZ66000014				
B3	KZ66000015				
C1	KZ66000016				
C2	KZ66000017				
C3	KZ66000018				
C4	KZ66000019				
D3	KZ66000020				
			Accessories		
			Durafill VS Shade Guide (16 Shades)	KZ66000046	
			1 PLT Applicator STEEL	KZ66019748	



WHITENING SOLUTIONS

Venus White® Max

Maximum whitening in a flash.

Remarkably efficient—accelerate patients’ whitening with results in an hour or less.

- Three or four 15-minute treatments with 38% hydrogen peroxide gel.
- Syringe-to-syringe mixing chemically activates the product, ensures a thorough mix and provides maximum freshness.
- Contains potassium nitrate which helps decrease sensitivity associated with bleaching.



VENUS WHITE® MAX	
Venus White Max In-Office Kit (2 x 0.96ml Hydrogen Peroxide, 2 x 0.24ml Activator, 1 x 1.2ml Gingival Barrier, 5 x 22ga Flock Tips, 2 x 20ga Mini Tips)	KZ40005211
Venus White Max Light Cure Gingival Barrier (4 x 1.2ml Syringes, 8 Black Mini Tips)	KZ40005202

Venus White® Pro

Custom whitening for exceptional impact.

High-performance whitening for custom trays take a patient’s smile up to seven shades whiter.*

- Choose from 16% or 22% carbamide peroxide strengths.
- Contains potassium nitrate which helps decrease sensitivity associated with bleaching.
- Mint flavored for patient comfort.

VENUS WHITE® PRO	
16% Patient Kit (6 x 1.2ml Syringes, Tray Case, Shade Guide)	KZ40005166
22% Patient Kit (6 x 1.2ml Syringes, Tray Case, Shade Guide)	KZ40005167
16% Refill Kit (3 x 1.2ml Syringes)	KZ40005164
22% Refill Kit (3 x 1.2ml Syringes)	KZ40005165
Venus Comfort Gel (4 x 1.2ml Syringes, 4 Tips)	KZ40005240

*Results may vary.

VENUS WHITE® ULTRA (HYDROGEN PEROXIDE GEL)	
11.2% Venus White Ultra (7 upper and 7 lower pre-filled trays)	KZ40005210

WHITENING SOLUTIONS

Venus White®
Whitening Toothpaste
Bright, healthy smiles maintained.

Venus White® Whitening Toothpaste helps maintain a bright smile while providing a fresh clean taste and anti-cavity protection.

- Venus White® Whitening Toothpaste provides anti cavity protection while gently removing surface stains.
- The low-abrasion formula provides maximum enamel and dental protection and leaves your mouth feeling clean.
- Contains natural and artificial flavors and ingredients.



VENUS WHITE® WHITENING TOOTHPASTE	
Venus White Teeth Whitening Toothpaste - 4oz	KZ66051675
Venus White Teeth Whitening Toothpaste - 1oz	KZ66051676

Venus White® Teeth Whitening Touch Up Brush
Whenever, wherever whitening.

The Venus White® Teeth Whitening Touch Up Brush is a convenient, reusable pen filled with mint-flavored whitening gel, and is an easy way to help maintain a white smile.

- Portable design makes it ideal for on-the-go touch ups.
- Use the mint-flavored, 9% hydrogen peroxide formula up to twice a day for two weeks.
- Great option to help maintain the results achieved from Venus White Ultra, Pro or Max.



VENUS WHITE® TEETH WHITENING TOUCH UP BRUSH	
Venus White Teeth Whitening Touch Up Brush	KZ66051761



RESTORATIVE SOLUTIONS

High bond strength for reliable and consistent results

iBOND® Universal

8th generation - Our definition of Universal bonding is convenience.



iBOND Universal is the universal, light-curing adhesive for any bonding technique and all indications. Being compatible with all dental materials, it is the all-rounder for reliable bond strength.

Benefits

- Easy, convenient handling.
 - short curing time (10s)
 - efficient, controlled dispensing and application via drop control bottle design
 - can be stored at room temperature
 - usable directly out of bottle (no mixing)
- No change of personal habits
 - can be used in the preferred bonding technique (self etch, selective enamel etch or total etch (etch & rinse))
- Low technique sensitivity
 - one bottle (except iBond Ceramic Primer for glass/silicate ceramics)
- Instant bonding success/reliable bond strength
 - fast evaporation of water
 - homogeneous layer
- All dental materials
 - for all indications, type of restorations and materials (for glass/silicate ceramics: iBond Ceramic Bond need to be used additionally)
- compatible to light-, dual- and self-curing resin cements and core build up materials without activator
- Allows intraoral repair
 - minimally invasive treatment to restore chippings, fractures or localised marginal caries lesions instead of expensive replacement
- All indications
- All dental materials
- Multi-technique capabilities
- Easy handling
- Drop control
- Proven instant bonding success

Indications

- Bonding of direct restorations for all cavity classes (Black using light-curing, dual-curing or self-curing methacrylate-based composites/ compomers).
- Bonding of light-curing, dual-curing or self-curing core build-up materials.
- Sealing of hypersensitive tooth areas
- Sealing of cavities prior to amalgam restorations
- Bonding of fissure sealants.
- Sealing of cavities and core preparations prior to temporary cementation of indirect restorations (according to the immediate dentine sealing technique).
- Cementation of indirect restorations with light-curing, dual-curing or self-curing adhesive resin cements.
- Intraoral repair of composite and compomer restorations, porcelain fused to metal, metal and all-ceramic restorations.



iBOND UNIVERSAL		
Product	Contents	Art. code
iBOND universal bottle refill	1 x 4ml bottle; Pictorial card	KZ66061411
iBOND Universal bottle value pack	3 x 4ml bottles; Pictorial card	KZ66061412
iBOND Universal single dose value pack	100 Single doses; 50 application tips green (slim); 50 application tips red (regular) Pictorial card	KZ66061419
iBOND ceramic primer	1 x 4ml bottle	KZ66061416

RESTORATIVE SOLUTIONS

iBOND® Total Etch
Consistent reliability.

High bond strength for reliable and consistent results

iBOND Total Etch is a etch and rinse 2-step adhesive that attains high bond strength to both enamel and dentin while ensuring optimal marginal sealing. It is also ideal for desensitizing hypersensitive areas.

- The next generation of Gluma Comfort Bond & Desensitizer.
- Clinically proven to provide reliable and consistent bond strength.
- Optimized with nano-filler for improved intrinsic strength.



iBOND TOTAL ETCH

iBOND Total Etch Bottle Refill (1 x 4ml Bottle, 50 Application Tips, IFU and Pictorial Card)	KZ66040094
iBOND Total Etch Single Dose Refill (50 Single Doses, 50 Application Tips, IFU, and Pictorial Card)	KZ66040093
iBOND Total Etch Bottle Assortment (1 x 4ml Bottle, 50 Application Tips, Mixing Well, 2 x 2.5ml iBOND Etch 35 Gel Syringes, 25 Cannulas)	KZ66039866
iBOND Etch 35 Fluid (1 x 15ml Bottle iBOND Etch 35 Fluid)	KZ66039868
iBOND Etch 35 Gel (2 x 2.5ml iBOND Etch 35 Syringes, and 25 Cannulas)	KZ66039862

Gluma® Desensitizer PowerGel
Spot on desensitization.

With the control of a gel and confidence of Gluma, your patients will appreciate the fast, effective, and long-term results for a more comfortable prophy.

- Gel consistency for greater control and accuracy during placement.
- No-drip formula, so it stays where you place it.
- Green color indication for easy visual placement and rinse-clean capability.



GLUMA® DESENSITIZER POWERGEL

Gluma Desensitizer PowerGel Syringe Kit (4 x 1gm)	KZ66043451
Gluma Desensitizer PowerGel Brush Cannulas	KZ66044334



RESTORATIVE SOLUTIONS

Translux® Wave®

For safe, sustainable polymerisation.

Translux® Wave is a LED curing light for the polymerisation of dental materials such as adhesives and composites.

It is developed to consistently produce high-quality fillings with an optimum polymerisation depth in the 440 – 480 nm wavelength range.

- Pen-style body: The light-weight and ergonomic pen-style design enables an effortless handling even for small hands. Furthermore, the advanced angle of the 360° rotatable light guide allows a better access to the molar regions.
- One-button operation: One convenient power/program button to choose easily between 10 and 20 seconds. LED lights indicate the program change. provides freedom of movement
- Cordless design: The cordless design gives you maximum freedom of movement.
- Reliable performance: Translux Wave comes with a strong and long life Li-ion battery. It requires very little time to recharge and the battery can be exchanged quickly and easily without tools. The built-in radiometer ensures the performance you need for the perfect restoration.



TRANSLUX WAVE	
Order Information	Art. code
Hand piece, 8mm light guide, glare shield, charging unit, power supply	KZ66055013
Accessories	Art. code
Light guide, 8mm, 70°	KZ66057365
Rechargeable battery	KZ66057363
Glare shield, orange	KZ66016455
Replacement contacts for charger	KZ66057382
Hand Piece incl. battery	KZ66057383
Charging unit grey	KZ66057384
Power supply	KZ66057468

TECHNICAL DATA	
Power Density	> 1,200 mW/cm² (Cure Rite)
Wavelength	440 – 480nm (narrow)
Dimensions charging unit	L140 x W60 x H60mm
Dimensions hand piece	L210 mm, Ø 23mm
Weight	Charging unit 437g, Hand piece: 145g
Power Supply (charging unit)	100 to 240 V(AC), 50/60 Hz, 5 W, 12 – 18VA
Power Supply (hand piece)	Rechargeable Lithium-ion battery 3.6 V/2000mAh
Light Source	High-luminosity mono-die LED
Light guide	Ø 8 mm, 70°

RESTORATIVE SOLUTIONS

Translux® 2Wave®

Made to perform.

Translux® 2Wave® the light-weight cordless LED curing light with a 360° rotatable light beam, offers maximum freedom of movement and technical novelties. We developed its features to perform for your benefit:

- Maximum freedom of movement and minimized hand fatigue thanks to the light weight of only 150 g and the cordless pen-style design.
- The light guide, having a high angle of 70° and being 360° rotatable, makes the molar region easily accessible.
- Due to the LED display and simple controls handling becomes intuitive.
- The scratch-resistant surface is quickly and easily disinfected.
- The built-in radiometer enables you to check the performance of the curing light.
- The Translux 2Wave has strong, long-lasting Li-ion batteries that need only very little time for recharging.
- The Translux 2Wave features a steady high-intensity output and a homogenous light beam that are ideal for curing all common light-curing composites and adhesives.
- The broad spectral output with two wavelength ranges enables you to cover new and innovative photoinitiators as well as common camphor-quinone initiators.



TRANSLUX 2WAVE	
Order Information	Art. code
Hand piece, 8mm light guide, glare shield, charging unit, power supply	KZ66055013
Accessories	Art. code
Light guide, 8mm, 70°	KZ66057365
Rechargeable battery	KZ66057363
Glare shield, orange	KZ66016455
Replacement contacts for charger	KZ66057366
Hand Piece incl. battery	KZ66057371
Charging unit grey	KZ66057372
Power supply	KZ66057468

TECHNICAL DATA	
Power Density	> 1,400 mW/cm² (Cure Rite)
Wavelength	385 – 510nm (narrow)
Dimensions charging unit	L140 x W60 x H60mm
Dimensions hand piece	L210 mm, Ø 23mm
Weight	Charging unit 440g, Hand piece: 150g
Power Supply (charging unit)	100 to 240 V(AC), 50/60 Hz, 12 – 18VA
Power Supply (hand piece)	Rechargeable Lithium-ion battery 3.6 V/2000mAh
Light Source	High-luminosity quad-die LED
Light guide	Ø 8 mm, 70°

RESTORATIVE SOLUTIONS

Ivory® Rubber Dam Clamps

The most popular stainless steel clamps.



IVORY® MOST POPULAR

Complete Rubber Dam Kit (Punch, Forceps, Rubber Dam, Frame, 8 Clamp Starter Kit, and Template)	KZ50057966
Rubber Dam Clamp Starter Kit (8 Clamps: 0, 2, 2A, 7, 8A, 9, 14, 14A)	KZ50057968

ANTERIOR CLAMPS

Stainless Steel Clamp No. 6	KZ50057382
Stainless Steel Clamp No. 9	KZ50057384

TIGER CLAMPS

Stainless Steel Clamp No. 1T	KZ50057841
Stainless Steel Clamp No. 2T	KZ50057843
Stainless Steel Clamp No. 2AT	KZ50057847
Stainless Steel Clamp No. 9T	KZ50057855
Stainless Steel Clamp No. 14T	KZ50057849
Stainless Steel Clamp No. 14DT	KZ50057857
Stainless Steel Clamp No. 56T	KZ50057853

- Available in a multitude of sizes and shapes
- Die-cut, heat-treated, tempered and individually hand-set to ensure high performance
- Also available in wingless

WINGED CLAMPS

Stainless Steel Clamp No. 00	KZ50057302
Stainless Steel Clamp No. 0	KZ50057304
Stainless Steel Clamp No. 1	KZ50057306
Stainless Steel Clamp No. 2	KZ50057310
Stainless Steel Clamp No. 2A	KZ50057312
Stainless Steel Clamp No. 3	KZ50057322
Stainless Steel Clamp No. 4	KZ50057324
Stainless Steel Clamp No. 7	KZ50057328
Stainless Steel Clamp No. 8	KZ50057334
Stainless Steel Clamp No. 8A	KZ50057336
Stainless Steel Clamp No. 12A	KZ50057348
Stainless Steel Clamp No. 13A	KZ50057352
Stainless Steel Clamp No. 14	KZ50057354
Stainless Steel Clamp No. 14A	KZ50057356
Stainless Steel Clamp No. 56	KZ50057374
Stainless Steel Clamp No. 1A	KZ50057308
Stainless Steel Clamp No. 5	KZ50057326
Stainless Steel Clamp No. 7A	KZ50057330
Stainless Steel Clamp No. 10	KZ50057338
Stainless Steel Clamp No. 11	KZ50057342
Stainless Steel Clamp No. 26	KZ50057370
Stainless Steel Clamp No. 27	KZ50057372



IVORY® RUBBER DAM

Premium Rubber Dam	
Blue 6" x 6"—36 Sheets	KZ40000266
Green 6" x 6"—36 Sheets	KZ40000281

ACCESSORIES

Rubber Dam Forceps	KZ50057220
Rubber Dam Punch	KZ50057225

Rubber Dam Products	
Adult Frame (6 in)	KZ50057230

Rubber Dam Products	
Template 6" x 6"	KZ50057235

Flexitime®



IMPRESSION SOLUTIONS

Your Clinical Problem Solver*

A-Silicone Impression Material Flexitime® VPS

Precise detail for accurate fit every time.

Flexitime combines over 15 years of market experience with the outstanding performance and features of a precision impression material.

- Reliable impressions even in adverse conditions.
- Superior hydrophilicity on polyether level.
- Exceptional toughness prevents tearing and distortion.

FLEXITIME® – VPS IMPRESSION MATERIAL

Fast&Scan Easy Putty (1 x 600ml)	KZ66045754
Fast&Scan Heavy Tray (2 x 50ml Cartridges with Mix Tips)	KZ66060792
Fast&Scan Light Flow (2 x 50ml Cartridges with Mix Tips)	KZ66044871
Fast&Scan Medium Flow (2 x 50ml Cartridges with Mix Tips)	KZ66044881
Easy Putty (1 x 600ml)	KZ50034802
Heavy Tray (2 x 50ml Cartridges with Mix Tips)	KZ50034804
Monophase (2 x 50ml Cartridges with Mix Tips)	KZ50034811
Medium Flow (2 x 50ml Cartridges with Mix Tips)	KZ40005041
Correct Flow (2 x 50ml Cartridges with Mix Tips)	KZ50034806
Light Flow (2 x 50ml Cartridges with Mix Tips)	KZ40005040

ACCESSORIES

Mix Tips 4.2mm Yellow (48 Count)	KZ66001988
Mix Tips 6.5mm Green (48 Count)	KZ66001989
Intraoral Tips, Yellow (96 Count Small)	KZ66000783
Universal Silicone Adhesive (10ml)	KZ65607758
Dispensing Gun 2 (2:1/1:1)	KZ66001970

FLEXITIME 15 YEARS



IMPRESSION SOLUTIONS

Flexitime® Dynamix® speed Automatic Mixing Machine

The perfect mix. Every time.

Dynamix speed offers all of the advantages of automatic mixing, most importantly it saves the practice time.

- Void free mixing ensures precision and accuracy.
- Working time: 60-seconds to 2:30 minutes; Total set time: 3:30 to 5 minutes.



FLEXITIME® DYNAMIX® SPEED

Dynamix Fast&Scan Heavy Tray (2 x 380ml)	KZ66045755
Dynamix ProScan Monophase refill (2 x 380ml)	KZ66045801
Dynamix Heavy Tray Refill (2 x 380ml)	KZ40003991
Dynamix Monophase Refill (2 x 380ml)	KZ40003990
Dynamix Putty Refill (2 x 380ml)	KZ40003992

ACCESSORIES

Dynamix Mixing Tips (1 x 50 Mixing Tips)	KZ40004353
Dynamix Fixation Ring (1 x 2 count)	KZ66047385



FLEXITIME® DYNAMIX® SPEED AUTOMATIC MIXING MACHINE

Dynamix speed EU Standard Automatic Mixing Machine	KZ66045779
--	------------

IMPRESSION SOLUTIONS

C-Silicone Impression Material – Optosil® / Xantopren®

Our state-of-the-art C-silicones.

The Optosil/Xantopren product assortment is proven in practice since many years and has set the standard for silicone precision impression early on. It enables precise reproduction of detail and shows no sensitivity either towards latex gloves or retraction solutions. The Optosil/Xantopren system is available in different versions tailor-made for your professional demands.

Optosil P Plus and Xantopren are reliable for daily use in the dental practice. As a matched system of classical putty and low viscosity materials, they are user friendly and enable precise reproduction of detail. As such they are our classics for daily use in the dental practice.

C-SILICONE IMPRESSION MATERIAL

Optosil Comfort 1 x 900ml	KZ50034202
Optosil Liquid Combi	KZ66037715
Optosil P Plus 1 x 900ml	KZ66046475
Xantopren Comfort Light 2 x 50ml	KZ50034102
Xantopren Comfort Medium 2 x 50ml	KZ50034105
Xantopren Function Combi Pack	KZ66037442
Xantopren H Green 1 x 140ml	KZ66046476
Xantopren L Blue 1 x 140ml	KZ66046448
Xantopren M Mucosa 1 x 140ml	KZ66046447
Xantopren Mixing Cup	KZ65278223
Xantopren Spatula	KZ65283677
Xantopren VI Plus 1 x 140ml	KZ66046449
Periphery Wax 60 Sticks	KZ50092189
Copr Wax Bite	KZ50092178
Activator Universal Plus Liquid 25ml	KZ66065537
Activator Universal Plus Paste 60ml	KZ66037443



IMPRESSION SOLUTIONS

Xantasil® Alginate substitute

Impressively convenient.

Xantasil covers all anatomical impression indications for which traditional alginates would typically be used. It can be used as a substitute for traditional alginates for taking impressions of temporaries, orthodontic models, splints and impressions of opposing jaws.



XANTASIL® ALGinate SUBSTITUTE

Xantasil Cartridge Fast Set 6 x 50ml	KZ66052102
Xantasil Cartridge Fast Set 2 x 380ml	KZ66052101

Flexitime® Bite

The foundation for accurate results.

Flexitime Bite delivers exceptional accuracy for standard bite registrations and superior powder free scannability and optimization.

- Sets fast in just 30 seconds.
- Extremely high shore D hardness of 40.
- Whipped cream like consistency ensures exceptional accuracy.



FLEXTIME® BITE

Flexitime Bite Single
(2 x 50ml Cartridges and 6 Mix Tips) **KZ50066038**

**Scannable
without
powder!**

Memosil® 2 / Memoreg® 2

Transparent VPS bite registration material.

High degree of dimensional stability with accurate and consistent results.

- Ideal material for fabricating standard and light-cured provisional laminates, crowns and bridges.



MEMOSIL 2

Memosil 2 (2 x 50ml Cartridges and 12 Mix Tips)	KZ50035018
Memoreg 2 (2 x 50ml Cartridges and 12 Mix Tips)	KZ50035517

ACCESSORIES

Dispensing Gun 2 (2:1/1:1)	KZ66001970
Mix Tips 6.5mm Green (48 Count)	KZ66001989
Universal Silicone Adhesive (10ml)	KZ65607758

IMPRESSION SOLUTIONS

ALGINOPLAST®

The Firm Alginate Offering detailed Reproduction with high stability.


It shows a high level of accuracy, elasticity and tensile strength, and is also extremely compatible with plaster. The dust free and clean mixing is easy and comfortable.

We use high quality seaweed from preferred cultivation areas as a basis for our alginate impression material. These alginates last longer without preservatives. Since they do not contain any chemical disinfectants, they are skin and gingiva friendly. Their mild peppermint taste also makes the impression more pleasant for the patient. Alginoplast is best suited for documentation, analysis, planning, opposing bite and working model.

Alginoplast fulfils your requirements of modern materials for anatomic impression

- Its stable consistency, offering detailed reproduction, is ideal for dentists who prefer a stable alginate with high stability.
- The fast setting variant saves time and is practical.
- The normal-setting variant is suitable for automatic mixing without time pressure during processing.

ALGINOPLAST®	
Alginoplast Nh 1 x 500g	KZ66046478
Alginoplast Sh 1 x 500g	KZ66046446



XANTALGIN®

The Fine-Paste Alginate For Sensitive Patients.

It shows a high level of accuracy and is particularly suitable for sensitive patients due to its soft starting consistency. The dust free and clean mixing is easy and comfortable.

We use high quality seaweed from preferred cultivation areas as a basis for our alginate impression material. These alginates last longer without preservatives. Since they do not contain any chemical disinfectants, they are skin and gingiva friendly. Their mild peppermint taste also makes the impression more pleasant for the patient. Xantalgin is best suited for documentation, analysis, planning, opposing bite and working model.


Xantalgin Select

- The creamy consistency is particularly good for sensitive patients.
- Thanks to its high elasticity and tensile strength, impressions can be removed easily, even if there was undercutting.
- The short setting time is pleasant and saves time.

Xantalgin Crono

- Xantalgin Crono shows each phase of the setting reaction through a colour change.
- Thanks to its optimised reaction mechanism, the impressions can be poured out up to 5 days later.

XANTALGIN®	
Xantalgin Crono 1 x 500g	KZ66046450
Xantalgin Select 1 x 500g	KZ66046477
Xantalgin Mixing Cup	KZ65234064
Xantalgin Spatula	KZ65234072



ACCESSORIES

VibraKleen E2® Ultrasonic Cleaning Tablets

The natural power of enzymes.

- Powerful protease enzymes quickly remove and breakdown hard to clean deposits like blood, bone, and food particles that chemical detergents can leave behind.
- Powerful, yet gentle on your equipment.
- Biodegradable and pH neutral.
- Non-toxic, non-corrosive formula contains no chlorine or harsh alkaline to damage your expensive stainless steel equipment of temporaries, orthodontic models, splints and impressions of opposing jaws.




VIBRAKLEEN E2®	
VibraKleen E2 / 80 Tabs	KZ50037101

Venus® Supra

Long-lasting, high gloss polishers.

Use **Venus® Supra** for polishing all types of composites - in just two steps to achieve long-lasting, high gloss polishing effects for naturally beautiful composite restorations. Additionally, Venus Supra polishers achieve extremely low surface roughness, reducing plaque accumulation and pigmentation over time.

- Silicone polishers highly filled with microfine diamond powder.
- Available in 5 different shapes: Cup, Disk, Flame, Large Cup, Large Flame.
- Autoclavable and Latex-free.



VENUS® SUPRA	
Venus Supra Kit (10 polishers: each of all five shapes as pre- and high gloss polishers)	KZ66039110
Refill-Cup (10 polishers: One shape as pre- and high gloss polishers (each 5x))	KZ66039161
Refill-Disk (10 polishers: One shape as pre- and high gloss polishers (each 5x))	KZ66039162
Refill-Flame (10 polishers: One shape as pre- and high gloss polishers (each 5x))	KZ66039163
Refill-Large Cup (10 polishers: One shape as pre- and high gloss polishers (each 5x))	KZ66039164
Refill-Large Flame (10 polishers: One shape as pre- and high gloss polishers (each 5x))	KZ66039165

Sopira® Citoject® Syringes

Anaesthesia injections with pinpoint precision.

Less discomfort for your patients and a more relaxing working atmosphere for you. SOPIRA Citoject makes intraligamentary injections more comfortable – an outstanding innovative design available in stainless steel or light alloy.

- Balanced weight distribution for stable handling and precision.
- The integrated plastic sleeve offers added safety.
- Unparalleled control for injection.



SOPIRA CITOJECT	
SOPIRA Citoject Barr Str 2.2ml/inch LA	KZ66002396
SOPIRA Citoject Plastic Sleeve 1.8ml 1pCE	KZ66005419
Sopira Citoject Plastic Sleeve 2.2ml x 1	KZ66006866
Sopira Citoject Set 2.2ml/Inch Flat LA	KZ66002395



KULZER

MITSUI CHEMICALS GROUP

ADAM DENTAL
supplies

1300 449 262

sales@adamdental.com.au
17 Pickering Rd, Mulgrave VIC 3170, Australia
www.adamdental.com.au