

PROMOTION
OCT-DEC 2022

# **PRIMELINE**





□ AeroAcoustic blade



Glass rod fiber optics



Powerful & Precise Power Head - 25W



Reverse flow protection



Automatic pressure regulation 1



Incredibly Quiet Only 55 dB noise



Refined ceramic bearings



Scratch resistance Iconic Prime body



Coupling Connection Kavo

25 Watt HP21KL HP21NL HP21WL

HP21SL

Power Head

Small Head 18 Watt HP22KL

HP22NL HP22WL - N/A -

\$1,683.00



Now with 2 Years FairCare Warranty<sup>2</sup>

Unconditional warranty except physical damage on Prime Line High-speeds

#### **V** COUPLINGS

NSK

W&H

Sirona









MK-dent & Kavo with LED & water

\$385.00



 QC5016W \$385.00



Sirona connection

with LED & water

 QC6016KW \$495.00

QC6016NT

with fibre optic

 QC6016SW \$495.00

Non-optic

 ✓ Non-optic QC4014NT \$296.00

❷ Non-optic QC4014W \$296.00

Non-optic -NA-

 QC4014K \$296.00

1) Only available for Kavo coupling 2) FairCare warranty applicable on Prime Line high speeds sold after 1st Oct 2018

- Buy 3 Get 15% OFF
- Buy 4 Get 20% OFF







2 Years FairCare Warranty









- % 1:5 Red •
- O Max. 200,000 rpm
- O Internal water
- For FG-burs (Ø 1.6 mm)
- LP15L (F/Optic) RRP: \$1,826.00

Buy 3 - \$1,552.10 ea

Buy 4 - \$1,460.80 ea

- % 1:1 Blue •
- O Internal water
- For CA-burs (Ø 2.35 mm)
- LP11L (F/Optic) RRP: \$1,287.00

Buy 3 - \$1,093.95 ea

Buy 4 - \$1,029.60 ea

- 🗞 5.4:1 Green •
- Max. 7,407 rpm ○ Internal water
- For CA-burs
- (Ø 2.35 mm)
- LP54L (F/Optic) RRP: \$1,320.00

Buy 3 - \$1,122.00 ea

Buy4 - \$1,056.00 ea

- % 1:1 Blue •
- C Max. 40,000 rpm
- O Internal water
- For HP-burs (Ø 2.35 mm)
- LP01L (F/Optic)RRP: \$1,045.00

Buy 3 - \$888.25 ea

Buy 4 - \$836.00 ea

AIR MOTORS 🕏





- LED
- AM1006 \$990.00



- **AM1004** \$748.00





• **AM0004** \$ 440.00

### NEW

### KAVO COMPATIBLE HEADS V



- Kavo 68LH compatible
- Blue (1:1)
- LCH11 \$ 440.00



- Kavo 67LH compatible
- Green (2:1)
- LCH21 \$ 440.00



- Kavo 31LR compatible
- Green (2:1)
- LCH21P \$ 440.00



3 x Eco Line HS with F/O for any coupling connection

Savings from \$426



62 dB Operating Noise



AeroPerformance blade



Refined ceramic bearings

Jpgrade to 24W



Glass rod fiber optics



Smooth Eco body





3 x Eco Line HS F/O for any connection **PLUS** 

3 x Blue (1:1) Basic Line Contra-Angles LB11/11





Fibre Optic

5 x Eco Line HS F/O for any connection **PLUS** 

3 x Blue (1:1) Basic Line Contra-Angles LB11/11



Now with 1 Year FairCare Warranty Unconditional warranty except physical damage



### Order codes:

Coupling	Power Head	Small Head	Power Head	Small Head
Connection	24W (F/Optic)	16W (F/Optic)	24W (N/Optic)	16W (N/Optic)
Kavo	HE21KL	HE22KL	HE21K	HE22K
NSK	HE21NL	HE22NL	HE21N	HE22N
W&H	HE21WL	HE22WL	HE21W	HE22W
Sirona	HE21SL	– N/A –	– N/A –	– N/A –

• \$890.00 (F/Optic)

• \$649.00 (N/Optic) Buy 3 - Get 10% OFF

### **♥** BASIC LINE CONTRA - ANGLE HEADS (NSK-EC Compatible) -



- For CA burs (Ø 2.35mm)
- LBH11 \$231.00



- For CA burs (Ø 2.35mm)
- LBH11C \$165.00



- For Prophy Snap-on cups
- LBH11P \$198.00

<sup>1)</sup> FairCare warranty (except physical damage) applicable on Eco Line sold after 1st Jan 2019 2) Eco Line Power Head 24 Watt

- Buy 3 Get 15% OFF
- Buy 4 Get 20% OFF
   Mix & Match



### FairCare Warranty

I year unconditional warranty except physical damage



- 🗞 1:5 Red 🛭
- $\bigcirc$  Internal water
- For FG-burs (Ø 1.6 mm)
- LE15L (F/Optic) RRP: \$1,485.00
- Buy 3 \$1,262.25 ea

### Buy 4 - \$1,188.00 ea

- LE15 (Non-optic) RRP: \$1,287.00
- Buy 3 \$1,093.95 ea

Buy 4 - \$1,029.60 ea



- 🗞 1:1 Blue •
- O Internal water
- For CA-burs (Ø 2.35 mm)
- LE11L (F/Optic) RRP: \$990.00
- Buy 3 \$ 841.50 ea

### Buy 4 - \$ 792.00 ea

- LE11 (Non-optic)
   RRP:\$770.00
- Buy 3 \$ 654.50 ea

Buy 4 - \$ 616.00 ea



- § 5.4:1 Green ●
- O Internal water
- Screw-in & Snap-on Prophy Cups
- LE54LP (F/Optic) RRP: \$1,210.00
- Buy 3 \$ 1,028.50 ea

### Buy 4 - \$ 968.00 ea

- LE54P (Non-optic) RRP: \$935.00
- Buy 4 \$ 748.00 ea



% 1:1 Blue ●

○ Internal water

For HP-burs

(Ø 2.35 mm)

• LE01L (F/Optic)

- Buy 4 \$ 748.00 ea

   LE01 (Non-optic)
- RRP: \$748.00
- Buy 3 \$ 635.80 ea

Buy 4 - \$ 598.40 ea

BASIC LINE CONTRA - ANGLES V



- % 1:1 Blue ●
- O Push button
- LB11/11
- Buy 1 \$418.00
- Buy 3 \$339.00 ea



- 🗞 1:1 Blue •
- O Latch Type
- LB11/11C
- Buy 1 \$418.00
- Buy 3 \$339.00 ea



- % 4:1 Green •
- O Snap-on Cups
- LB41/11P
- Buy 1 \$ 462.00
- Buy 3 \$396.00 ea



- % 1:1 Blue ●
- LB02
- Buy 1 \$ 418.00
- Buy 3 \$339.00 ea



### Buv 3 - Get 1 FREE





MK-deut PROPHY LINE



EMS Non-optic Type EN-041

 SC21E Buy 1 - \$440.00 Buy 3 - Get 1 FREE

**EMS Piezon Non-optic** Type EN-061

 SC21E2 Buy 1 - \$539.00 Buv 3 - Get 1 FREE



**EMS Piezon LED** Type EN-060A

 SC21EL Buy 1 - \$715.00 **Buv 3 - Get 1 FREE** 



KaVo Piezo LED Type 1.007.3995

 SC21KL\* NEW Buy 1 - \$715.00

Buy 3 - Get 1 FREE \*Only for use with EMS compatible scaler tips

Acteon-Satelec Non-optic Type F 12 281

SC21S

Buy 1 - \$ 440.00

Buy 3 - Get 1 FREE



Acteon-Satelec LED Type F 12 609

 SC21SL Buy 1 - \$715.00

Buy 3 - Get 1 FREE

**♥** SCALER TIPS-



ST01SD

\$88.00

SATELEC 1

Buy 3 - Get 1 FREE

- High-quality Diamond-Like-Carbon coating
- Extremely hard and wear-resistant surface

### • SCALING DLC COATED Only available for Satelec fitting

ST02SD SATELEC 1S \$88.00

ST03SD **EMSA** \$88.00

ST04SD SATALEC 2 \$132.00

PERIO DLC COATED Only available for Satelec fitting

ST06SD SATELEC 10Z \$132.00

ST10SD SATALEC 10P \$154.00

ENDO DLC

ST08SD SATALEC S07 \$154.00

SCALING STAINLESS STEEL Only available for Satelec fitting

ST02S SATELEC 1S \$63.80

ST03S **EMSA** \$63.80

ST04S SATELEC 2 \$88.00

PERIO STAINLESS STEEL

ST07E/ST07S **EMS DPL3** SATELEC 10X \$110.00

ST09E/ST09S **EMS HPL** 

\$110.00

ST10E/ST10S SATALEC 10P

\$110.00

ST11E/ST11S EMS RT1

\$110.00

EMS and Acteon-Satelec are registered trademarks. MK-dent GmbH has no economic connection with the aforementioned companies.

# **Universal Torque Wrench**

### TOZQ CONTZOL®

Precision instrument with a Torque Adjustment System for tightening of Prosthetic components at specific Torque, for any implant system.



# Automatic Declutching when Selected Torque is reached

- Compliance with Protocols for screw tightening
- Limits the risk of the screw fracturing or loosening



### Universal connection of Prosthetic mandrels

- 7 different Torque levels (from 10 to 35 N.cm)
- Can be used with ALL Implant Systems



# Micro-Head and Ergonomic Handle

 Easy access in restricted area with Excellent Visibility



### Easy and Quick to use

 Reduced motion compared with using a ratchet wrench while tightening




Lightweight ONLY 135 gms

Order code: 15501

### 3 Easy Steps









www.anthogyr.com





 Ergonomically shaped, compact and light-weight air-driven polishing device 360° rotatable outer handle due to unique RotaGrip® for simple and fatique-proof operation



- Prophy Powder (300 gms) \$60.50 each
  - Lemon PP1012

### Additional Prophy tip



SP3011 \$374.00

### Coupling connection order codes:

Kavo - PR1011K, NSK - PR1011N, W&H - PR1011W, Sirona - PR1011S

### V LEDs -



- Quality illumination, No Glare
- 5500K Natural Daylight Colour
- Clear Vision & Sharper Contrast
- Longer Durability & Economical



\$88.00

BU8012

BU8012KM

- Kavo & MK-dent motors

BU8012WH - W&H handpieces

BU8012WLT - W&H airmotors LT25 BU8012WEA - W&H motors 40LT, 50 LT & A51LT

BU8012SG - Sirona coupling and motors

BU8012N

NSK coupling

BU8012NT

NSK motors TiM40 & Ti205L

Kavo & MK-dent coupling

BU8012NNL - NSK NL motors

BU8012NM - NSK coupling KCL,SCL,PTL LED

#### ♥ PREMIUM LUBRICATION OIL



Fully Synthetic

For any brand of handpieces and motors

LU1011 (500 ml) \$68.20



Fully Synthetic

For Kayo Quattrocare 2 lubrication system

LU1022 (500 ml) \$ 68.20



For W&H Assistina **lubrication** system

 LU1031 (500 ml) \$ 88.00

MANUFACTURER

MX-dont MK-dent GmbH

Marie-Curie-Str. 2 22941 Bargteheide

**GERMANY** 

www.mk-dent.com info@mk-dent.com YOUR SPECIALIST DEALER